






Choose with confidence

Motiva Round

	MINI 			DEMI 			FULL 			CORSEÉ 		
 (cm)	$\tilde{D}_{\sim B}$ (cm)	$\tilde{D}_{\sim C}$ (cm)	V (cc)	$\tilde{D}_{\sim B}$ (cm)	$\tilde{D}_{\sim C}$ (cm)	V (cc)	$\tilde{D}_{\sim B}$ (cm)	$\tilde{D}_{\sim C}$ (cm)	V (cc)	$\tilde{D}_{\sim B}$ (cm)	$\tilde{D}_{\sim C}$ (cm)	V (cc)
8.5	2.2	5.0	105	3.1	5.7	135	3.5	6.0	145	4.0	6.5	180
9	2.3	5.3	125	3.3	6.0	155	3.7	6.3	175	4.2	6.8	210
9.5	2.4	5.5	140	3.4	6.3	180	3.9	6.7	205	4.5	7.2	240
9.75	2.4	5.6	150	3.4	6.4	190	4.0	6.9	220	4.6	7.4	260
10	2.5	5.8	160	3.5	6.5	205	4.1	7.1	235	4.8	7.7	280
10.25	2.5	5.9	170	3.5	6.6	215	4.2	7.2	255	4.9	7.9	300
10.5	2.6	6.1	185	3.6	6.8	230	4.3	7.4	275	5.1	8.2	325
10.75	2.6	6.2	205	3.7	7.0	245	4.4	7.6	295	5.2	8.3	350
11	2.7	6.4	220	3.8	7.1	265	4.5	7.7	315	5.4	8.6	380
11.25	2.7	6.5	230	3.8	7.2	285	4.6	7.9	335	5.5	8.7	410
11.5	2.8	6.6	245	3.9	7.4	300	4.7	8.1	355	5.7	9.0	440
11.75	2.8	6.7	260	3.9	7.5	320	4.8	8.2	375	5.8	9.2	475
12	2.9	6.9	275	4.0	7.7	340	4.9	8.4	400	6.0	9.5	510
12.25	2.9	7.0	290	4.0	7.8	360	5.0	8.6	425	6.1	9.6	550
12.5	3.0	7.2	310	4.1	7.9	380	5.1	8.8	450	6.3	9.9	590
13	3.1	7.5	360	4.3	8.3	425	5.3	9.1	500	6.6	10.3	650
13.5	3.2	7.7	400	4.4	8.5	475	5.5	9.5	550	6.9	10.8	725
14	3.3	8.0	430	4.5	8.8	525	5.7	9.8	625	7.2	11.2	825
14.5	3.4	8.3	475	4.6	9.1	575	5.9	10.2	700	7.5	11.7	925
15	3.5	8.6	525	4.8	9.4	625	6.1	10.5	775	7.8	12.1	1050

A=Base B=Projection C=Arc Length V=Volume

SilkSurface® PLUS

Mini	Demi	Full	Corsé	ProgressiveGel PLUS™	BluSeal®	Single Pack	TwinPack™	Sizers	Q Inside	2Y Extended Warranty*
✓	✓	✓	✓	✓	✓	✓	✓	✓	✗	✓

SilkSurface® PLUS Q

Mini	Demi	Full	Corsé	ProgressiveGel PLUS™	BluSeal®	Single Pack	TwinPack™	Sizers	Q Inside	5Y Extended Warranty*
✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

VelvetSurface™ PLUS Q

Mini	Demi	Full	Corsé	ProgressiveGel PLUS™	BluSeal®	Single Pack	TwinPack™	Sizers	Q Inside	5Y Extended Warranty*
✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

*Optional



www.motivaimplants.com